



HPLC columns


Connected chromatography solutions

Low-flow columns and accessories

Introduction

Low-flow chromatography is ideal when detailed sample information is required from small sample volumes, such as proteomic, metabolomic, and intact protein analysis. The Thermo Scientific range of nano-, capillary-, and micro-flow columns offer excellent sensitivity and resolution in easy-to-use formats.


- Thermo Scientific™ μ PAC™ Neo HPLC Columns
- Thermo Scientific™ EASY-Spray™ HPLC Columns
- Thermo Scientific™ Double nanoViper™ HPLC Columns






Video: Low-flow HPLC columns connectivity


Contents

μ PAC Neo HPLC columns	5
<hr/>	
EASY-Spray HPLC columns	8
Bottom-up proteomics	9
Top-down proteomics	10
<hr/>	
Double nanoViper HPLC columns	12
Bottom-up proteomics	13
Top-down proteomics	14
<hr/>	



Column selection guide

	Pillar array column format	Packed bed column format	
	Thermo Scientific™ μ PAC™ Neo HPLC Columns	Thermo Scientific™ EASY-Spray™ HPLC Columns	Thermo Scientific™ Double nanoViper™ HPLC Columns
Technology			
Benefits	<p>Ultimate separation</p> <ul style="list-style-type: none"> • Excellent retention time stability • A unique combination of performance and reliability to get the highest sample coverage every time • Separate emitters • Compatible with all low-flow U/HPLC instruments 	<p>Ease-of-use</p> <ul style="list-style-type: none"> • Click-and-spray connect with Thermo Scientific™ EASY-Spray™ Source • Thermo Scientific™ nanoViper™ connections • Integrated column and emitter • Integrated temperature control • For use with Thermo Scientific™ mass spectrometry systems 	<p>Analytical flexibility</p> <ul style="list-style-type: none"> • Universal Thermo Scientific™ nanoViper™ Fingertight Fittings for column inlet and outlet • Simple zero-dead-volume (ZDV) connections • Separate emitters • Compatible with all low-flow U/HPLC instruments
Application areas/chemistries	<p>Deliver excellent column-to-column reproducibility with flow rate flexibility. Ideally suited for proteomic analyses of HPLC separations up to 450 bar.</p> <ul style="list-style-type: none"> • 50 cm column: 15–60 min gradient time • 110 cm column: 90–150 min gradient time • 50 cm Low Load single cell analysis: 15–60 min gradient time • High Throughput: <15 min gradient time, for up to 180 samples/day 	<p>Bottom-up proteomic applications</p> <p>The Thermo Scientific™ PepMap™ Neo UHPLC Columns are a recent addition to our portfolio. PepMap Neo columns are packed to higher pressure, which provides 1500 bar pressure rating, improved column-to-column consistency, and increased efficiency.</p> <p>Top- and middle-down proteomic applications</p> <p>The Thermo Scientific™ MAbPac™ Capillary Reversed-Phase HPLC Column is best suited for the characterization of intact proteins in top- and middle-down proteomic applications where sample amount is limited.</p>	



μPAC Neo HPLC columns



The μPAC Neo columns are specifically suited for bottom-up proteomic applications where separation performance is critical to the success of the analysis. Our μPAC Neo HPLC columns offer highest resolution and peak capacities for complex biological samples. The unique μ-pillar backbone improves column-to-column reproducibility and robustness, providing more confidence in analytical results.

Additional reading

Links	Type	Description
	Reference guide	Chromatography consumables reference guide for low-flow LC-MS proteomic research
	Flyer	Enabling high sensitivity LC-MS analysis for bottom-up and top-down proteomic research
	Learn more thermofisher.com/lowflowHPLCcolumns	

Choose a μPAC Neo HPLC column when:

- Highest resolution and peak capacities is required
- Your samples span a wide concentration range
- Highest LC-MS sensitivity is needed
- You want to speed up your runtimes
- LC-MS robustness is needed
- You want an increased column lifetime
- You prefer working at much lower back pressures than with packed bed columns
- It is important to compare results from experiments spanning over time or geographical location

What makes μPAC Neo HPLC columns special?

The unique separation path provides:

- μ-pillar stationary backbone, micromachined in a silicon wafer
- Flow path designed for highest analyte concentration during elution
- Extra high-resolution separations, using up to 110 cm column lengths
- Low back pressure separations, improving column and emitter robustness
- Perfect match with single cell proteomics sample amounts



μPAC Neo HPLC columns Continued

Velocity Label-free Quantitation (LFQ) Data Independent Acquisition (DIA) Platform




Figure 1. Graphical schematic of HR-DIA workflow for label-free quantitation of two- and three-proteome mixtures. The different components of the workflow are depicted on the top. The main goal of the setup is the quantitative performance at high sample throughput while delivering robust and reproducible results to make it a perfect fit for large scale clinical and biomarker discovery studies.




Figure 2. HR-DIA Workflow delivers confident proteome coverage utilizing next generation library-free analysis approaches. Bar graph comparison of protein group (human and *E. coli*) and precursor (total) numbers identified in 12 runs of two-proteome mix by use of three different software packages. Data analysis has been done by library-free analysis. All protein group results are filtered for 1% experiment-wide FDR.

μPAC Neo HPLC columns Continued



Ordering information


Description	Pillar dimensions (μm)	Interpillar distance (μm)	Column length (cm)	Flowrate range (nL/min)	Cat. no
Thermo Scientific 50 cm μPAC Neo Column	2.5	1.25	50	100-750	COL-NANO050NEOB
Thermo Scientific 110 cm μPAC Neo Column	2.5	1.25	110	100-750	COL-NANO110NEOB
Thermo Scientific 50 cm μPAC Neo Low Load column	2.5	1.25	50	100-750	COL-LOLO050NEOB
Thermo Scientific μPAC Neo High Throughput Column	3.0	2	5.5	250-2500	COL-CAPHTNEOB

Ordering information

Description	Pillar dimensions (μm)	Interpillar distance (μm)	Column length (cm)	Cat. no
Thermo Scientific™ μPAC™ Neo Trapping Column	5	2.5	1	COL-TRPLOLONEOB2

Ordering information

Description	Pillar dimensions (μm)	Details	For use with	Cat. no
Thermo Scientific™ EASY-Spray™ Nano Emitters	10	Bullet type without transfer line	EASY-Spray ion-source	ES993



EASY-Spray HPLC columns



Ensure robust nano- and capillary-flow LC-MS analysis using Thermo Scientific EASY-Spray HPLC Columns. The integrated column/emitter design eliminates dead volume and is temperature-controlled for maximum reliability and performance. Rigorously tested to ensure maximum quality, these columns deliver maximum simplicity and ease-of-use. The capillary-flow HPLC columns provide sensitive protein, peptide, and monoclonal antibody (MAb) separation. They give proteomic researchers more than ever before: more throughput, more sensitivity, more separation power, and more ease of use.

Additional reading

Links	Type	Description
	Reference guide	Chromatography consumables reference guide for low-flow LC-MS proteomic research
	Flyer	Enabling high sensitivity LC-MS analysis for bottom-up and top-down proteomic research
	Learn more thermofisher.com/lowflowHPLCcolumns	

Choose an EASY-Spray column when:

- You want simple connections with an EASY-Spray source. This is ideal for novice and experienced users
- Sample amount is limited
- Analytical UHPLC does not provide sufficient sensitivity
- Workflow simplicity is key
- High sensitivity is required to identify proteins and peptides at low expression levels
- Analyses are done in a targeted and untargeted way for screening and verification

What makes an EASY-Spray column special?

Unique design provides uncompromised performance in an ease-of-use format for nano and capillary LC-MS analysis

Features for optimum data quality:

- Simple connection to the LC and Thermo Scientific MS instruments
- Precision machined and positioned glass emitters
- Integrated nanoViper zero-dead-volume (ZDV) unions
- Integrated temperature control



Video:

Thermo Scientific EASY-Spray
150 mm LC columns





PepMap Neo HPLC columns

The Thermo Scientific™ EASY-Spray™ PepMap™ Neo UHPLC Columns are perfect for bottom-up proteomic applications. Packed at higher pressure and rated to 1500 bar, they provide consistent column-to-column performance, long column lifetime, and excellent efficiency. These benefits are true at any pressure.




Figure 3. The 60% reduction in total analysis time allows increasing the sample throughput moving from the nano- to the capillary-flow LC-MS method



Ordering information for bottom-up proteomic applications

Description	Length (mm)	Column ID (μm)	Cat. no
EASY-Spray PepMap Neo UHPLC Columns	150	75	ES75150PN
	500	75	ES75500PN
	750	75	ES75750PN





MABPac Capillary Reversed Phase HPLC Column

The Thermo Scientific™ MABPac™ Capillary Reversed Phase capillary column is best suited for the characterization of intact proteins in top-down proteomic, clinical, and anti-doping applications where sample amount is limited or sensitivity is crucial.




Figure 4. Calculation of site occupancy of N306 in Fab glycosylated mAb



Ordering information for top-down proteomic applications

Description	Length (mm)	Column ID (μm)	Cat. no
EASY-Spray HPLC Column	150	150	ES907



EASY-Spray HPLC columns Continued



EASY-Spray accessories

For the best performance from your EASY-Spray column consider investing in these accessories.




Ordering information

Description	Union type	Particle size (μm)	Column ID (μm)	Media bed length (mm)	Trap length (mm)	Cat. no
Thermo Scientific™ PepMap™ Neo Trap Cartridge	N/A	5	300	5	N/A	174500
	Nut/sleeve	5	200	20	150	164213
	Nut/sleeve	5	100	20	150	164199
Thermo Scientific™ Acclaim™ PepMap™ 100 C18 HPLC Column, Trap Column	Double nanoViper	5	100	20	150	164750
	Double nanoViper	3	75	20	150	164535
	Double nanoViper	3	75	20	70	164946

Ordering information

Description	For use with	Cat. no
Thermo Scientific™ PepMap™ Neo Trap Cartridge Holder, PEEK Tubing, and nanoViper™ Fittings	Low-flow PepMap columns	174502



Double nanoViper columns



The Thermo Scientific™ Viper™ and Thermo Scientific™ nanoViper™ Fingertight Fitting Systems provide tool-free connections designed to be used for the entire fluidic pathway in LC systems to improve chromatographic results.

Virtually without any dead-volume, Viper and nanoViper fittings combine usability with high performance. Viper and nanoViper connections can be used on all standard LC modules, valves, and columns quickly, independent of different connection geometries and system backpressures. Dedicated capillary kits for standard LC system configurations and application-specific setups enable high qualitative and reproducible results for all flow rates and pressure ranges.

Additional reading

Links	Type	Description
	Reference guide	Chromatography consumables reference guide for low-flow LC-MS proteomic research
	Flyer	Enabling high sensitivity LC-MS analysis for bottom-up and top-down proteomic research
	Product specifications	Viper and nanoViper Fingertight Fitting Systems
	Learn more thermofisher.com/lowflowHPLCcolumns	

Choose these columns when:


- Maximum flexibility is required
- Changing the emitter and column independently is important



What makes these columns special?


These stand-alone nano-, capillary-, and micro-flow columns are:

- Designed with single nanoViper and double nanoViper fingertight fittings for trouble-free connection
- For robust separation in proteomics research, drug discovery, and high-throughput proteomics laboratories



Video:

Discover a better LC connection



Double nanoViper PepMap Neo UHPLC columns

Separate challenging peptide mapping samples with Thermo Scientific™ Double nanoViper™ PepMap™ Neo UHPLC Columns. These columns feature easy connectivity, high reproducibility, and excellent separations. Our Neo columns are packed to higher pressure and provide 1500 bar pressure capability, improved column-to-column consistency, and increased efficiency. The column media is manufactured and selected to exacting standards and packed at high pressure, resulting in enhanced peak symmetry, resolution, and column-to-column reproducibility that allows you to obtain greater sample coverage and sample insights.




Figure 5. Reproducible identification and quantification of HeLa peptides and proteins over 4 EASY-Spray PepMap Neo columns while using the Vanquish Neo UHPLC system coupled with the Orbitrap Exploris 480 mass spectrometer



Ordering information for bottom-up proteomic applications

Format	Length (mm)	Column ID (µm)	Cat. no
Double nanoViper PepMap Neo UHPLC Columns	150	75	DENV75150PN
	500	75	DENV75500PN
	750	75	DENV75750PN

Double nanoViper columns Continued

Top-down proteomics



MABPac Capillary Reversed Phase HPLC Column

The Thermo Scientific MABPac Capillary Reversed Phase column is best suited for the characterization of intact proteins in top-down proteomic, clinical, and anti-doping applications where sample amount is limited or sensitivity is crucial.




Figure 6. Calculation of site occupancy of N306 in Fab glycosylated mAb



Ordering information for top-down proteomic applications

Format	Length (mm)	Column ID (μm)	Cat. no.
MABPac Capillary Reversed Phase HPLC Column	150	150	164947



Double nanoViper columns Continued



LC-MS connection accessories and emitters

These emitters, nanoViper tubing kits, and unions offer easy connection from your LC system to an EASY-Spray source.



Ordering information

Description	For use with	Part number
Viper and nanoViper Fingertight Fittings Accessories		6040.2304
nanoViper Fingertight Fittings, 20 µm x 550 mm	Double nanoViper columns	6041.5260
EASY-Spray Nano Emitter, 10 µm		ES993
EASY-Spray Capillary Emitter, 15 µm		ES994

Traps and accessories

For the best performance from your double nanoViper column consider investing in these nanotraps.



Ordering information

Description	Union type	Particle size (µm)	Column ID (µm)	Media bed length (mm)	Trap length (mm)	Cat. no
Thermo Scientific™ PepMap™ Neo Trap Cartridge	N/A	5	300	5	N/A	174500
	Nut/sleeve	5	200	20	150	164213
	Nut/sleeve	5	100	20	150	164199
Thermo Scientific™ Acclaim™ PepMap™ 100 C18 HPLC Column, Trap Column	Double nanoViper	5	100	20	150	164750
	Double nanoViper	3	75	20	150	164535
	Double nanoViper	3	75	20	70	164946



Expect reproducible results with sample prep, columns and vials



Don't see what you need? We would be happy to discuss your specific requirements. Please contact your local sales representative for custom orders.



Learn more at thermofisher.com/lowflowHPLCcolumns